TEST NAME/ PSYCHE CODE	COLLECTION GUIDELINES	REFERENCE VALUES	ROOM TEMP.	REFRIGERATED	FROZEN	Clinical Utility
SEDIMENTATION RATE (ESR)	1 EDTA LAVENDER TOP. (1 ml	MALE: 0-20 mm/HOUR	4 hours	12 hours		The erythrocyte sedimentation rate (ESR) is an easy, inexpensive, nonspecific test that has been used for
CODE: ESR	MINIMUM)	FEMALE: 0-30 mm/HOUR				many years to help detect conditions associated with acute and chronic inflammation, including infections,
CPT: 85651	SCHEDULE: DAILY					cancers, and autoimmune diseases.
SEMEN ANALYSIS  CODE: SEMEN	PERFORMED 7am-2pm, M-F. WILL NOT BE ACCEPTED AFTER 2 PM. DELIVER TO LAB WITHIN 30 MINUTES OF COLLECTION, KEEP AT BODY TEMPERATURE.	INTERPRETIVE REPORT				A semen analysis is used to determine whether a man might be infertile
CPT: 89320	TEIVIPERATURE.					
SEMEN (COUNT ONLY)						This test is used to evaluate a man's fertility. It can also be used to confirm the success of a vasectomy.
CODE: SPRM	PERFORMED 7am-2pm M-F. <u>WILL NOT</u> <u>BE ACCEPTED AFTER 2 PM</u> . DELIVER TO THE LAB WITHIN 30 MINUTES, KEEP AT BODY TEMPERATURE.	> THAN 40 MILLION.				
SENSITIVITY, ANTIBIOTIC (MIC)	AUTOMATICALLY ORDERED WITH			1 week		Sensitivites may be added on as requested. Example: yeast, fungal on wound or urine cultures.
CODE: N/A	PATHOGEN GROWTH.					
CPT: 87186	SCHEDULE: DAILY					
SODIUM, SERUM	1.0 ml SERUM	135-148 mmol/L		48 hours		Blood sodium testing is used to detect abnormal concentrations of sodium
CODE: NA	1.0 IIII SEKOIVI	155-146 ПППОІ/С				
CPT: 84295	SCHEDULE: DAILY					
SEX HORMONE BINDING GLOBULIN	1.0 ML SERUM	MALE: 20-50YR 13.3-89.5nmol/L	8 hours	7 days	2 months	The test for sex hormone binding globulin (SHBG) is primarily ordered in conjunction with other tests to
CODE: SHB		FEMALE: 20-46YR 18.2-135.5nmol/L				evaluate the status of a person's androgen levels - the male hormones.
		POST MENOPAUSE: 47-91YR 16.8- 125.2nmol/L				
CPT: 84270	SCHEDULE: DAILY					
SODIUM, URINE (RANDOM AND 24 HOUR)	12 ml ALIQUOT OF 24 HOUR OR	RANDOM: NO RANGE ESTABLISHED 24 HOUR:50-260 meq/24 HOURS	2 hours	24 hours		Sodium results are used to diagnose and monitor conditions involving electrolyte inbalance
CODE: RUNA (RANDOM)	RANDOM URINE.					
CODE: 24UNA						
CPT: 84300	SCHEDULE: DAILY					